

Development and testing of an inhalable dust sampling system for 20 l/min

PEROSH 2025 – 6th Research Conference Manchester

PEROSH, Manchester

10.09.2025

Dörr, Marcel

Table of content

1. Aim
2. Specifications of standards
3. Aerosol chamber experiments
4. Wind channel experiments
5. Measurement uncertainty
6. Summary / Outlook



1. Aim

- Validation of 4 inhalable dust prototypes for personal aerosol sampler for 20 l/min.
- Generate reproduceable results
 - dust chamber and wind tunnel needs to be adjusted



UFO1



UFO2



Beetle



Cone


2. Specifications of standards: EN 13205

Variable	Terms of use	Measurement conditions
Aerodynamic particle diameter	Inhalable dust: 10 μm – 100 μm	≥ 9 different particle diameter, 1 needs to be $\geq 90\mu\text{m}$
Wind velocity	Indoor and outdoor workplaces	$\leq 0,1$ m/s and 1m/s
Wind direction	Personal sampling	Continuous rotation or stepwise ≥ 4 directions
Variability of sampler test samples	Test group as large as possible	≥ 6 Sampler
Relative standard deviation of dust conc.	At the measurement plane	$\sigma \leq 10 \%$

Additionally (EN482)

- Target concentration should be 0.1x to 2.0x the OELV, 10 mg/m^3 inhalable dust


2. Specifications of standards: EN 481




3. Aerosol chamber



3. Aerosol chamber – Sampling efficiency at 0,1 m/s, 20 mg/m³





3. Aerosol chamber – Sampling efficiency at 0,1 m/s, 1 mg/m³




4. Wind tunnel



4. Wind tunnel – Measurement setup, 1 m/s



4. Wind tunnel – Sampler positioning setup



5. Measurement uncertainty - Overview

Source	Measurement uncertainty	
Dust mass losses due to filtration	F1200 (10 µm)	2,5%
	F240 (90 µm)	15,5%
Inaccuracy of the scales	absolute	0,01 mg
Sample adhesion in the sampler	GSP 3,5: ~1%	GSP 10: ~2%
	UFO1: 1,9%	UFO2: t.b.d.
	Beetle: 2,04%	Cone: 2,69%
Relative standard deviation of dust conc.	Aerosol chamber	3,97%
	Wind tunnel	9,01%
Whole measurement uncertainty	t.b.d.	

6. Summary / Outlook

- Aerosol chamber
 - 5 to 9 dusts measured
 - High conc.: sampling efficiency did not lower with the particle size
 - Low conc.: sampling efficiency follows EN481
 - For choosing 1 prototype more information necessary
- Wind tunnel
 - Determination of the dust for which the selected setup complies with the standard ($\sigma < 10\%$)
 - Due to higher sedimentation of coarser dusts
 - Modify the nozzle device, install a traverse or move it manually

Thank you for

your attention

PEROSH, Manchester

10.09.2025

Dörr, Marcel